Reinigen von Laseroptiken

Cleaning Laser Optics

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The correct cleaning of laser optics plays a major role in ensuring their high quality. The cleanliness of a substrate is very important, particularly when handling high power levels. Otherwise, the laser beam could destroy the coating and with it the optic.

We are happy to share a wealth of experience in the cleaning of laser optics with our customers and have therefore put together a few notes on the correct handling of said optics.

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Tips on Cleaning

- **Dusting off Optics**
Larger dust particles or sand grains should always be cleaned up with a bellows or nitrogen. They should absolutely never be wiped off, as the substrate can become scratched in this way. Blowing a substrate with your mouth can cause miniature droplets of saliva to land on the surface.
- **Wearing Gloves When Handling Optics**
When handling optics gloves should be worn at all times. If fingerprints end up on the substrate, they should be removed immediately. The residue left by fingerprints contains aggressive acids that continue to dissolve the surface as time goes on. Even after just a short while, the residue can no longer be removed.
- **Metal Coatings Cannot be Cleaned**
Metal coatings should never be cleaned. These coatings are extremely soft. Even lightly wiping the optics can leave scratches in the coating.
- **Cleaning Wipes for Glasses Should not be Used**
The cleaning wipes that are used on glasses are not suitable for cleaning laser optics. They often contain defogging substances that could leave residues on the optic and harm the coating layers, thereby reducing the laser damage threshold.

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Cleaning Instructions

An acetone and isopropanol resistant cloth is used for cleaning. It is folded (1) to create a plane surface for cleaning and (2) so that the cloth can be held with self-closing tweezers. Please be sure that the metal part of the tweezers does not come into contact with the substrate to prevent the surface being cleaned from becoming scratched.

The cloth is lightly moistened with acetone or isopropanol. The optimal amount can be determined by experimenting. Too much solvent can leave streaks and too little makes the cloth become fuzzy too quickly.

The substrate can best be cleaned by starting in the middle and working in circular motions out to the edge. A constant, not too slight pressure should be applied with the cloth to the surface being cleaned.

The optic should be cleaned in transmission over a cold light source and then inspected.

As soon as the cloth begins to become fuzzy, it must be replaced by a new one.

Standard safety requirements should be heeded when handling cleaning agents. This includes, for example, wearing solvent resistant gloves, using a well-ventilated room, and heeding all hazard sheets.

For more information on cleaning, please feel free to contact us at any time.