


Datasheet

OE-300-IN-01

200 MHz Variable Gain Photoreceiver



The image shows model OE-300-IN-01-FC.


Features

- Adjustable transimpedance gain from  $10^2$  to  $10^8$  V/A
- Wide bandwidth up to 200 MHz
- InGaAs-PIN photodiode covering the 900 to 1700 nm wavelength range
- FC fiber optic input
- High dynamic input range up to 10 mW optical power
- Very low noise, NEP down to 47 fW/√Hz
- Switchable low pass filters for minimizing wideband noise
- Full manual and remote control capability

Applications

- All-purpose low-noise photoreceiver (O/E converter) for the MHz range
- Time resolved optical pulse and power measurements
- Laser intensity noise measurements (RIN)
- Optical front-end for oscilloscopes, spectrum analyzers, A/D converters and RF lock-in amplifiers

Block Diagram






BS-OE-300-R1

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY



DE-OE-300-IN-01\_R2/MG/JM/250CT2018

<b>Datasheet</b>		<b>OE-300-IN-01</b>
<b>200 MHz Variable Gain Photoreceiver</b>		
Available Versions	OE-300-IN-01-FC	FC fiber optic input
Related OE-300 Models	See separate datasheets for following models on <a href="http://www.femto.de">www.femto.de</a> :	
	OE-300-SI-10-FST	Si-PIN, 1 mm x 1 mm, 400 - 1000 nm 1.035"-40 threaded flange
	OE-300-SI-10-FS	Si-PIN, 1 mm x 1 mm, 400 - 1000 nm 25 mm dia. unthreaded flange
	OE-300-SI-30-FST	Si-PIN, $\varnothing$ 3 mm, 320 - 1000 nm 1.035"-40 threaded flange
	OE-300-SI-30-FS	Si-PIN, $\varnothing$ 3 mm, 320 - 1000 nm 25 mm dia. unthreaded flange
	OE-300-IN-03-FST	InGaAs-PIN, $\varnothing$ 300 $\mu$ m, 800 - 1700 nm 1.035"-40 threaded flange
	OE-300-IN-03-FS	InGaAs-PIN, $\varnothing$ 300 $\mu$ m, 800 - 1700 nm 25 mm dia. unthreaded flange
	OE-300-S	customized versions available on request
Available Accessories	<p><b>PRA-PAP</b>  post adapter plate, easy to mount on FEMTO photoreceiver series OE, FWPR, HCA-S and LCA-S</p> <p><b>PS-15</b>  power supply, input: 100 - 240 VAC, output: <math>\pm</math>15 VDC, +400/-250 mA</p> <p><b>LUCI-10</b>  compact digital I/O interface for USB remote control, supports opto-isolation of amplifier signal path from PC USB port, 16 digital outputs, 3 opto-isolated digital inputs, bus-powered operation</p>	
<b>SOPHISTICATED TOOLS FOR SIGNAL RECOVERY</b>		<b>F E M T O</b>

**Datasheet**

**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

**Specifications**

Test conditions  $V_s = \pm 15\text{ V}$ ,  $T_a = 25\text{ }^\circ\text{C}$ , system impedance = 50  $\Omega$

**Gain**

Transimpedance gain  $1 \times 10^2 \dots 1 \times 10^8\text{ V/A}$   
Gain accuracy  $\pm 1\%$

**Frequency Response**

Lower cut-off frequency DC/100 Hz, switchable  
Upper cut-off frequency up to 200 MHz (see table below),  
switchable to 1 MHz or 10 MHz

**Input**

Noise equivalent power (NEP) see table below  
Max. CW saturation power see table below

**Detector**

Detector InGaAs-PIN photodiode  
Active area Integrated ball lens, suitable for fibers up to  
50  $\mu\text{m}$  core diameter  
  
Spectral response 900 - 1700 nm  
Sensitivity R 0.95 A/W typ. @ 1550 nm  
Dark current 0.02 nA typ.

**Performance Depending  
on Gain Setting**

Gain setting (low noise) (V/A)	$10^2$	$10^3$	$10^4$	$10^5$	$10^6$	$10^7$
Upper cut-off frequency (-3 dB)	200 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz
NEP ( $\sqrt{\text{Hz}}$ , @ 1550 nm)	180 pW	22 pW	1.9 pW	390 fW	140 fW	50 fW
Measured at	20 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz
Integrated input noise (RMS)*	4.9 $\mu\text{W}$	380 nW	23 nW	3.3 nW	0.84 nW	71 pW
CW sat. power (@ 1550 nm)	10 mW	1.0 mW	100 $\mu\text{W}$	10 $\mu\text{W}$	1.0 $\mu\text{W}$	100 nW

Gain setting (high speed) (V/A)	$10^3$	$10^4$	$10^5$	$10^6$	$10^7$	$10^8$
Upper cut-off frequency (-3 dB)	175 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz
NEP ( $\sqrt{\text{Hz}}$ , @ 1550 nm)	132 pW	6.3 pW	1.4 pW	350 fW	113 fW	47 fW
Measured at	18 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz
Integrated input noise (RMS)*	3.0 $\mu\text{W}$	285 nW	21 nW	3.2 nW	0.84 nW	71 pW
CW sat. power (@ 1550 nm)	1.0 mW	100 $\mu\text{W}$	10 $\mu\text{W}$	1.0 $\mu\text{W}$	100 nW	10 nW

\* The integrated input noise is measured with a shaded input in the full bandwidth ("FBW") setting (referred to 1550 nm). The measurement bandwidth is 3 x the upper cut-off frequency at the specific gain setting; filter slope is a 1<sup>st</sup> order roll-off.

The input referred peak-peak noise can be calculated from the RMS noise as follows:

$$P_{\text{input noise peak-to-peak}} = P_{\text{input noise RMS}} \times 6$$

The output noise is given by:

$$U_{\text{Output noise RMS}} = P_{\text{input noise RMS}} \times \text{gain} \times R$$

$$U_{\text{Output noise peak-to-peak}} = U_{\text{Output noise RMS}} \times 6 = P_{\text{input noise RMS}} \times \text{gain} \times R \times 6$$

The integrated noise will be reduced considerably by setting the low pass filter to "1 MHz" or "10 MHz" instead of "FBW". This is especially useful for continuous wave (CW) measurements.

**Output**

Output voltage range  $\pm 1\text{ V}$  (@ 50  $\Omega$  load), for linear amplification  
Output impedance 50  $\Omega$  (designed for 50  $\Omega$  load)  
Slew rate 1000 V/ $\mu\text{s}$   
Max. output current  $\pm 40\text{ mA}$   
Output offset compensation adjustable by offset potentiometer and external control  
voltage, output offset compensation range min.  $\pm 100\text{ mV}$

**Ext. Offset Control**

Control voltage range  $\pm 10\text{ V}$   
Offset control input impedance 15 k $\Omega$

**SOPHISTICATED TOOLS FOR SIGNAL RECOVERY**



**Datasheet**

**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

Specifications (continued)

Indicator LED	Function	overload						
Digital Control	Control input voltage range	LOW bit: -0.8 ... +1.2 V, HIGH bit: +2.3 ... +12 V						
	Control input current	0 mA @ 0 V, 1.5 mA @ +5 V, 4.5 mA @ +12 V						
	Overload output	non active: <0.4 V @ 0 ... -1 mA active: typ. 5 ... 5.1 V @ 0 ... 2 mA						
Power Supply	Supply voltage	±15 V						
	Supply current	+110/-90 mA (depends on operating conditions, recommended power supply capability min ±200 mA)						
	Stabilized power supply output	±12 V, max. 20 mA, +5 V, max. 150 mA						
Case	Weight	320 g (0.74 lb.)						
	Material	AlMg4.5Mn, nickel-plated						
DC Monitor Output	Monitor output gain	<table border="1"> <thead> <tr> <th>Mode</th> <th>Monitor gain</th> </tr> </thead> <tbody> <tr> <td>Low noise</td> <td>Gain setting divided by -1</td> </tr> <tr> <td>High speed</td> <td>Gain setting divided by -10</td> </tr> </tbody> </table>	Mode	Monitor gain	Low noise	Gain setting divided by -1	High speed	Gain setting divided by -10
Mode	Monitor gain							
Low noise	Gain setting divided by -1							
High speed	Gain setting divided by -10							
	Monitor output polarity	inverting						
	Monitor output voltage range	±1 V (@ ≥1 MΩ load)						
	Monitor output bandwidth	DC ... 1 kHz						
	Monitor output impedance	1 kΩ (designed for ≥1 MΩ load)						
Temperature Range	Storage temperature	-40 ... +80 °C						
	Operating temperature	0 ... +60 °C						
Absolute Maximum Ratings	Max. CW power (averaged)	12 mW						
	Digital control input voltage	-5 V/+16 V relative to digital ground DGND (pin 9)						
	Analog control input voltage	±15 V relative to analog ground AGND (pin 3)						
	Power supply voltage	±20 V						

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY




**Datasheet**

**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

**Connectors**

Input FC fiber optic receptacle  
 Output BNC jack (female)  
 Power supply Lemo® series 1S, 3-pin fixed socket  
 (mating plug type: FFA.1S.303.CLAC52)  
 Pin 1: +15 V  
 Pin 2: -15 V  
 Pin 3: GND



Control port Sub-D 25-pin, female, qual. class 2  
 Pin 1: +12 V (stabilized power supply output)  
 Pin 2: -12 V (stabilized power supply output)  
 Pin 3: AGND (analog ground for pins 1 - 8)  
 Pin 4: +5 V (stabilized power supply output)  
 Pin 5: digital output: overload (referred to pin 3)  
 Pin 6: DC Monitor output  
 Pin 7: NC (= not connected)  
 Pin 8: output offset control voltage input  
 Pin 9: DGND (ground for digital control pins 10 - 16)  
 Pin 10: digital control input: gain, LSB  
 Pin 11: digital control input: gain  
 Pin 12: digital control input: gain, MSB  
 Pin 13: digital control input: AC/DC  
 Pin 14: digital control input: high speed / low noise  
 Pin 15: upper cut-off frequency limit 10 MHz  
 Pin 16: upper cut-off frequency limit 1 MHz  
 Pin 17 - 25: NC (= not connected)

**Scope of Delivery**

OE-300-IN-01-FC, Lemo® 3-pin connector, datasheet, transport package

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY



**Datasheet**


**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

**Remote Control Operation**

<b>General</b>	Remote control input bits are opto-isolated and connected by a logical OR function to the local switch settings. For remote control set the corresponding local switches to "Remote", "DC", "L" (low noise mode) and "FBW", and select the desired setting via a bit code at the corresponding digital inputs. Mixed operation, e.g. local AC/DC setting and remote controlled gain setting, is also possible.				
<b>Gain setting</b>	Low noise Gain (V/A) Pin 14=LOW	High speed Gain (V/A) Pin 14=HIGH	Pin 12 MSB	Pin 11	Pin 10 LSB
	10 <sup>2</sup>	10 <sup>3</sup>	LOW	LOW	LOW
	10 <sup>3</sup>	10 <sup>4</sup>	LOW	LOW	HIGH
	10 <sup>4</sup>	10 <sup>5</sup>	LOW	HIGH	LOW
	10 <sup>5</sup>	10 <sup>6</sup>	LOW	HIGH	HIGH
	10 <sup>6</sup>	10 <sup>7</sup>	HIGH	LOW	LOW
	10 <sup>7</sup>	10 <sup>8</sup>	HIGH	LOW	HIGH
<b>AC/DC setting</b>	<u>Coupling</u>	<u>Pin 13</u>			
	DC	LOW			
	AC	HIGH			
<b>Low pass filter setting</b>	<u>Upper cut-off freq. limit</u>	<u>Pin 15</u>	<u>Pin 16</u>		
	full bandwidth	LOW	LOW		
	10 MHz	HIGH	LOW		
	1 MHz	LOW	HIGH		
<b>High speed / low noise setting</b>	<u>Mode</u>	<u>Pin 14</u>			
	low noise mode	LOW			
	high speed mode	HIGH			

**Spectral Responsivity**



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY



**Datasheet**

**OE-300-IN-01**


**200 MHz Variable Gain Photoreceiver**

Typical Performance  
Characteristic


Frequency response

$$V_{\text{Supply}} = \pm 15 V_{\text{DC}}; R_{\text{Load}} = 50 \Omega$$


Gain setting: L10<sup>2</sup>, H10<sup>3</sup>




Gain setting: L10<sup>3</sup>, H10<sup>4</sup>




Gain setting: L10<sup>4</sup>, H10<sup>5</sup>




Gain setting: L10<sup>5</sup>, H10<sup>6</sup>



Gain setting: L10<sup>6</sup>, H10<sup>7</sup>



Gain setting: L10<sup>7</sup>, H10<sup>8</sup>



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY



**Datasheet**


**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**


Typical Performance  
Characteristic (continued)

Input noise equivalent power (NEP)


Gain setting:  $L10^2, H10^3$




Gain setting:  $L10^3, H10^4$




Gain setting:  $L10^4, H10^5$




Gain setting:  $L10^5, H10^6$



Gain setting:  $L10^6, H10^7$



Gain setting:  $L10^7, H10^8$



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY





**Datasheet**


**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

Typical Performance  
Characteristic (continued)


Signal pulse response

Gain setting L10<sup>2</sup>




Rise: 1.84 ns Fall: 1.90 ns

Gain setting H10<sup>3</sup>




Rise: 2.27 ns Fall: 2.32 ns

Gain setting L10<sup>3</sup>




Rise: 3.30 ns Fall: 3.41 ns

Gain setting H10<sup>4</sup>




Rise: 3.44 ns Fall: 3.52 ns

Gain setting L10<sup>4</sup>



Rise: 26.42 ns Fall: 26.49 ns

Gain setting H10<sup>5</sup>



Rise: 26.77 ns Fall: 27.01 ns

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY




**Datasheet**

**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

Typical Performance  
Characteristic (continued)



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY




**Datasheet**

**OE-300-IN-01**

**200 MHz Variable Gain Photoreceiver**

Dimensions

Fiber optic input OE-300-IN-01-FC:



Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany

**SOPHISTICATED TOOLS FOR SIGNAL RECOVERY**

