
A-CUBE-SERIES

Manual

1. DESCRIPTION OF THE APD-MODULE

1.1. General Information

A-CUBE is a highly sensitive APD-module which contains the following three main components: Avalanche photodiode, voltage source and transimpedance amplifier. As a plug-and-play module the A-CUBE features an integrated electrical circuit which can be powered by a voltage of ± 12 V.

Delivery contents:

- / A-CUBE module
- / Analogue adapter cable (Fig. 1)
- / Digital adapter cable (Fig. 2)

Optional delivery contents:

- / Coaxial cable (SMB-BNC)
- / Power supply unit »CUBE-PSU«

These instructions imply handling of the product with proper care and serve to improve the understanding of the working principle of the module. Included are the explanation of the usage and safety aspects related to the product.


Fig. 1: Adapter cable analogue
(for use with CUBE-PSU)


Fig. 2: Adapter cable digital
(for use with CUBE-PSU)

1.2. Technical Information

Two different types of A-CUBEs are available, containing either a Si-APD or an InGaAs-APD. Depending on the wavelength that is to be detected either one of these A-CUBE can be chosen:

APD-Type	Wavelength Range [nm]
Si	400 – 1100
InGaAs	1000 – 1700

The bandwidth is dependent on the gain. Therefore, several different types of A-CUBE are offered. Due to the differences in bandwidth the modules vary in terms of sensitivity. Further information can be found in the A-CUBE datasheet.

	1 MHz	3 MHz	10 MHz	25 MHz	100 MHz	180 MHz
Si-APD @ 905 nm [MV/W]	500	50	5	0.5	27	10
InGaAs-APD @ 1550 nm [MV/W]	94	9.4	0.94	0.094	5	1.9

2. USAGE OF THE A-CUBE

The A-CUBE is a plug-and-play module and is ready for use with a power supply unit and an oscilloscope.

The special housing of the module is designed to be suitable for optical experiments. Therefore, the A-CUBE can be mounted directly into a 30 mm cage system. After being put into the cage system the module can be locked in place using M2 set screws.

Using the additional drills (diameter 6 mm) on the front and the back side of the module more rods can be mounted, again using M2 set screws. Like this the A-CUBE can be integrated fully into the cage system.

Other mounting possibilities are available using the M6 screws at the bottom side of the module (Fig. 3).


Fig. 3: A-CUBE housing

2.1. Usage According to Regulations

No special measures for the protection of the module need to be taken.

The A-CUBE is only to be operated within the specifications mentioned in the datasheet.

The A-CUBE is not to be opened.

The customer simply needs to connect the adapter cable analogue and the power supply (with + 12 V, GND, -12 V) to the A-CUBE module. Alternatively, the customer can connect the digital adapter cable and the power supply to the A-CUBE and the PC. The power supply unit »CUBE-PSU« is recommended to be used with the module.

A coaxial cable (SMB-BNC) is to be connected to the A-CUBE (SMB) and a measuring device, for instance an oscilloscope (BNC).

2.2. Usage Not According to Regulations

Using supply voltages $V_1 > +15\text{ V}$ and $V_2 < -15\text{ V}$ can damage the electronics of the A-CUBE.

For damages that arise due to usage not according to regulations there is no product liability claim from LASER COMPONENTS Germany GmbH.

3. UNDERLYING GUIDELINES AND NORMS

EMV-Guideline: Guideline 2004/108/EG

RoHS 2 Guideline: Guideline 2011/65/EU

DIN EN 61326-1; VDE0843-20-1:2013-07

IEC 61326-1:2012;

German version: EN61326-1:2013

4. SAFETY INSTRUCTIONS

Data on the documents attached is binding.

Operating the A-CUBE is not dangerous for health and life.

5. SETUP AND FUNCTION OF A-CUBE MODULE

In the red and black housing there is the detector mounted in a TO-8 can with integrated amplifier and cover glass, the DC/DC-converter for the bias voltage of the APD and the logical circuit which is controlled by a microprocessor. At the back side of the housing there is an SMB-conconnector for the signal output, a 5-pin connector for the voltage supply connection (+ 12V, GND, -12V) and the connectors for the RxD and TxD interface. Moreover, there are two potentiometers where settings for bias voltage and RTC can be set manually (Fig. 4).

Control Element	Description
P1	Potentiometer to set U_{Bias}
P2	Potentiometer to set RTK
M9	Connector for supply voltage and RS232-communication
Signal Output (SMB)	Connector for signal output

Allocation of M9 Connector

Pin	Type	Designation
1	GND	Ground
2	+V _{op}	Operating voltage +12V/300mA
3	RxD	Recieve input (RS232)
4	-V _{op}	Operating voltage -12V/300mA
5	TxD	Transmit output (RS232)

By using the digital adapter cable without the CUBE-PSU, please note that you have to shorten the 2 GND Pins of Lumberg Connector KV 50 (Pin1 and Pin 2) to the same GND potential.

Allocation of M9 Connector Without CUBE-PSU

Pin	Type	Designation
1	GND	Ground
2	GND	Ground
3	n.c.	not connected
4	-V _{op}	Operating voltage -12V/300mA
5	+V _{op}	Operating voltage +12V/300mA


Fig. 4: Connectors of the A-CUBE


Fig. 5: Connector of Lumberg KV 50

6. STARTING-UP THE A-CUBE

The module is delivered in a protection packaging along with the data sheet and the manual. To start up the A-CUBE the protection packaging needs to be removed.

The A-CUBE can be connected to the power supply unit with the analogue adapter cable. In this case only the manual operation remains (see section 7.1).


The A-CUBE can be connected to the power supply unit and the PC with the digital adapter cable. In this case both the manual and the digital operation is possible (see section 7.1 and 7.2). As soon as the module is connected to the power supply, the module is active. By default the module is to be operated manually at first.

6.1. Digital Operation Setup - Setup of the Serial Interface in the Terminal Program:

Serial Port	
Baud rate	9600
Data	8 bit
Parity	none
Stop	1 bit
Flow control	none


Example HTerm

Example Tera Term

6.2. Setup for Receiving and Transmitting Data:

Terminal Setup	
Receive	CR
Transmit	CR


Example HTerm


Example Tera Term

6.3. Example for Command and Response of an A-CUBE

Transmitted Byte: 0100 1011 (0x4B => ASCII: K)

7. OPERATION OF THE A-CUBE

7.1. Manual Operation

The manual operation is active as soon as the module is powered up.

7.1.1. Setting the High Voltage

By turning the potentiometer P1 (see Fig. 4) the bias voltage of the APD can be set. By turning the potentiometer clockwise the bias voltage is increased, by turning it counterclockwise the bias voltage is decreased.

7.1.2. Setting the Resistor Temperature Compensation (RTC)

The RTC can be set by turning the potentiometer P2 (see Fig. 4). By turning the potentiometer clockwise the RTC is increased, by turning it counterclockwise the RTC is decreased. By turning P2 counterclockwise all the way the RTC is zero.

The RTC helps vary the voltage U_{Bias} depending on the temperature. This is necessary because the breakdown voltage of the APD varies with temperature as well.

If the temperature of the APD rises, the breakdown voltage of the APD rises as well and with it the operating point. If the temperature decreases the breakdown voltage of the APD decreases as well.

The factor with which the breakdown voltage varies depending on the temperature is the temperature coefficient (TC) of the APD. The TC needs to be determined individually for every APD and typically is 1V/°C for Si-APDs and 0.075V/°C for InGaAs-APDs.

7.2. Digital Operation

If the A-CUBE is connected to a PC and a power supply and a terminal programme is activated, the following appears in the terminal window (Fig. 6):


Fig. 6: Startup screen

The A-CUBE can be controlled via RS232 commands and provides the possibility of reading out the current data of the module. In order to do this, the manual operation mode needs to be switched to the digital operation mode via an RS232 command. This happens with the following commands:

RS232 Command	Function
r1	Switches from manual operation mode to digital operation mode
r2	Switches from digital operation mode to manual operation mode

7.2.1. Setting the Bias Voltage

RS232 Command	Function
uxxxx	Setting the bias voltage

xxxx stands for a value which is an element of the following set {0, ..., 4095}. The value which needs to be set corresponds to a value of the bias voltage. This value can be calculated using this formula:

$$xxxx \approx U_{\text{Bias}} [\text{V}] \cdot 8.75 [1/\text{V}];$$

It is important to note that xxxx needs to be an integer.

7.2.2. Setting the Compensation of the Temperature Coefficient (TC)

The compensation of the TC is set by using the following RS232 command:

RS232 Command	Function
txxxx	Setting the compensation of the TC

xxxx stands for a value which is an element of the following set {0, ..., 4095}. The value corresponding to the desired bias voltage can be calculated using:

$$U_{\text{Bias}} = U_{\text{Bias } 25^{\circ}\text{C}} + U_{\text{TK}} \cdot d_{\text{Temp}};$$

$$U_{\text{TK}} \cdot d_{\text{Temp}} = \text{xxxx} \cdot 0.0012 \text{ V} \cdot (T - 298 \text{ K});$$

7.2.3. Monitoring

Monitoring is possible in manual and digital mode.

The current values of bias voltage, temperature of the sensor and setpoints for both bias voltage and TC can be displayed using the following RS232 commands:

RS232 Command	Function
d3	continuous display of the following values (manual mode): / setpoint for U_{Bias} [mV] / setpoint for TC [mV] / current temperature [1/10 °C] / current U_{Bias} [1/10 V] continuous display of the following values (digital mode): / setpoint for U_{Bias} [digits] / setpoint for TK [digits] / current temperature [1/10 °C] / current U_{Bias} [1/10 V]
‚Esc‘	terminates the continuous display

7.2.4. Switching U_{Bias} on/off

The set bias voltage value can be switched on or off with the following RS232 command without losing the value or switching it back to 0:

RS232 Command	Function
m0	Setting the bias voltage to 0V
m1	Setting bias voltage to prior value

8. MAINTENANCE AND REPAIR

The A-CUBE module does not require any maintenance.

Contamination of the cover glass may lead to a reduction in sensitivity of the module. Therefore, please keep the cover glass clean. This can be done using a cloth suitable for cleaning glass. Please be careful when cleaning the glass not to leave scratches because this can reduce the sensitivity of the module as well.

9. DISPOSAL

The A-CUBE module was produced following the RoHS 2 guideline and can be disposed of with electrical waste.

10. TROUBLE SHOOTING

If no signal output can be detected please make sure if the beam hits the active area of the detector. Furthermore, please make sure via monitoring with the RS232 interface (see section 7.2.3.) that U_{Bias} is set correctly and is in contact with the APD.

11. WARRANTY


The warranty is limited to 12 months after delivery.

Repairs or changes with the module must be done by the manufacturer only. If this is not the case, the warranty of 12 months is invalid.

Further information on warranty is provided in the AGBs of LASER COMPONENTS Germany GmbH under the following link: <https://www.lasercomponents.com>

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